

# Examination Chair

## bon US-30

with powerful electromotive height adjustment  
and flip up armrest





## Ophthalmic examination chair US-30

The seat height of the US-30 is adjustable by electromotive control between 46 and 67 cm. It is operated conveniently via an existing examination unit or by a separate foot switch.


The examination chair is equipped with a folding armrest as standard. You can decide whether this should be placed to the left or right of the patient.

**CE** Like all **bon** units, also the US-30 examination chair is manufactured in Germany.

---

## Your new patient chair will be manufactured exactly to your wishes and needs.

You can choose between various useful features, such as a footrest in a fixed or folding version, a second swivelling armrest or a horizontal seat adjustment.



## Available in a variety of colours

The standard chair upholstery is made of an easy-care imitation leather. Your **bon** consultant will be pleased to show you our extensive selection of colour samples.

Alternatively, cloth or fine genuine leather is also available on request.