


**Amoils  
Epithelial  
Scrubber**



**Fast epithelial removal**

- Patient friendly
- **No knife or blade** used.
- Quick and efficient epithelial removal
- Less time for the corneal dehydration to occur.
- Produces easily seen "bull's eye" that aids in centration of laser application.
- Leaves smoother anterior stromal surface.
- Faster healing time, less epithelium removed.
- Particularly helpful for post RK/AK patients
- Preliminary studies suggest that **accuracy** of refractive results may be enhanced.

**Superior to Scrape**

- **2-5** seconds minimizes total procedure time
- Less patient anxiety by eliminating time consumed scraping
- Will not compromise stromal integrity due to top design.
- Eliminates nicks or scrapes to bowman's (Less possibility of keratocyte stimulation)
- Single use disposable tips, available as myopic and hyperopic

**Superior to Transepithelial**

- Uniform removal of the epithelium centrally and peripherally
- No need for subsequent scraping
- Less wear on laser optics
- Performed with patient positioned under the laser microscope
- Epithelial removal in re-treatment
- Essential in Hyperopic PRK and Corneal Crosslinking

**The Amoils Epithelial Scrubber Starter Kit**

includes:

- Power handle
- Neck
- Recharger
- 50 single use brush tips (25x myopic, 25x hyperopic)
- Fixation occluder
- Instructions with video